NOVISURGE

1 EHMS

Endoscopic Hernia Multifeed Stapler

EN	EHMS Endoscopic Hernia Multifeed Stapler
FR	Agrafeuse à alimentation multiple pour hernie par endoscopie EHMS
DE	EHMS Endoskopischer Multifeed Hernienhefter
ES	EHMS grapadora endoscópica para hernias de alimentación múltiple
PT	Agrafador Múltiplo Endoscópico para Hérnia EHMS
IT	Suturatrice endoscopica multipla per ernie EHMS
NL	EHMS endoscopische hernia multifeed stapler
DA	EHMS endoskopisk hernia (brok) stapler
SV	EHMS endoskopisk bräckhäftapparat med flermatning
FI	EHMS endoskoopinen tyräntöija
NO	EHMS endoskopisk multiavfyring-stiftmaskin for brokk
IS	EHMS holsjár fjölmötunarheftari fyrir kviðslit
TR	EHMS Endoskopik Fitik Çoklu Besleme Zimbasi
PL	Stapler endoskopowy wieloładunkowy do operacji przepuklin EHMS
CS	Endoskopický svorkovač hernií EHMS
RO	EHMS stapler endoscopic reîncărcabil pentru hernie
AR	دباسة EHMS مسبقة الماء لعلاج الفتق بالمنظار
HI	EHMS इंडोस्कोपिक हर्निया मल्टीफ़ीड स्टेपलर





Please consult the instructions for use carefully prior to use.

Failure to properly follow the instructions may lead to improper staple formation and, therefore, unsatisfactory stapling of tissue.

Important: This package insert is designed to provide instructions for use of the EHMS Endoscopic Hernia Multifeed Stapler. It is not a reference to surgical stapling techniques.

Instructions for use, Mode d'emploi, Gebrauchsanweisungen, Instrucciones de uso, Instruções de uso, Istruzioni per l'uso, Gebruiksaanwijzing, Brugsanvisning, Användarinstruktioner, Käyttöohjeet, Bruksanvisning, Leiðbeiningar um notkun, Kullanım talimatları, Instrukcja użytkowania, Návod k použití, Instrucțiuni de utilizare. تعليمات الاستخدام , उपयोग के लिए नर्दिदेश

Device Description

The EHMS Endoscopic Hernia Multifeed Stapler is a sterile, single patient use instrument designed for use with a 10/11 mm diameter surgical trocar. The device places rectangular titanium staples. Closed staple dimensions are approximately 5.3 mm x 3.5 mm.

Illustration and Nomenclature (Illustration 1)

1. Nose 2. Rotation Knob 3. Handle

Indications

The EHMS Endoscopic Hernia Multifeed Stapler has application for use in a variety of minimally invasive procedures for fixation of prosthetic material and approximation of tissue.

Contraindications

- This instrument is not intended for use in vascular or neural tissue or in solid organs.
- The staples in this instrument are not hemostatic and therefore should not be used to achieve hemostasis.

Instructions for Use

Verify compatibility of all instruments and accessories prior to using the instrument (refer to **Warnings and Precautions**).

1. Using sterile technique, remove the instrument from the package. To avoid damage, do not flip the instrument into the sterile field.
2. Introduce the instrument through a 10/11 mm diameter surgical trocar. (Illustration 2)
3. Under direct visualization, locate the site to be stapled.
4. Position the instrument over the desired location to be stapled using the rotation knob. For staple formation, squeeze the handle until reaching a full stop. (Illustration 3)
5. Release the handle completely to place the staple. (Illustration 4)
6. Inspect the site of staple placement to ensure complete staple formation and placement.
7. Repeat Steps 2-5 to place the next staple.

Warnings and Precautions

- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Dispose of all opened instruments whether used or unused.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

How Supplied and Stored

The EHMS Endoscopic Hernia Multifeed Stapler is supplied sterile and preloaded for single patient use and should be used prior to its expiry date identified on the box. Discard after use. Store product in cool, dry place.

Contact

Contact your local distributor or supplier to report problems or queries about the product.

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Changzhou Waston Medical Appliance Co., Ltd.
Advanced Surgical Solutions Group B.V.
Distributed by: Advanced Surgical Solutions Group B.V.